

## Use case table template

<b>Use Case Number:</b>	
<b>Use Case Name:</b>	
<b>Participating Actors:</b>	
<b>Goal:</b>	
<b>Trigger:</b>	
<b>Precondition:</b>	
<b>Postcondition:</b>	
<b>Basic Flow</b>	
<b>Exceptions</b>	
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	
<b>Open Issues:</b>	

<b>Use Case Number:</b>	1
<b>Use Case Name:</b>	<i>Start Skin Observer application</i>
<b>Participating Actors:</b>	main user (patient)
<b>Goal:</b>	initialize storage methods/check for previously stored data, open up the application for the user
<b>Trigger:</b>	this user clicks the app from the android app menu list
<b>Precondition:</b>	none
<b>Postcondition:</b>	The application is now at start up
<b>Basic Flow</b>	1/ User click application icon 2/ Application is launched
<b>Exceptions</b>	2.1/ application fails to launch display error message
<b>Qualities:</b>	System responds within 1s
<b>Functionality:</b>	<p>The application has:</p> <ul style="list-style-type: none"> <li>search/add/delete/view/edit</li> <li>---&gt; 1/ search: photo, skin condition, group</li> <li>---&gt; 2/ create: a photo, a skin condition(entry)</li> <li>---&gt; 3/ add: entries to group, photos to entries</li> <li>---&gt; 4/ delete: photos, entries, groups</li> <li>---&gt; 5/ view: photos, skin conditions, groups,</li> <li>---&gt; 6/ edit: photos, skin conditions, groups</li> </ul>
<b>Constraints:</b>	The UI looks aesthetically pleasing
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want to keep a visual log of any skin conditions I might . The purpose could be to gather a personal history of a medical condition as to provide evidence to a physician.
<b>Notes:</b>	<p>We'll add details into each function later</p> <p><b><u>-Tests for use case 1-</u></b></p> <ul style="list-style-type: none"> <li>1/ check if the application launches properly</li> <li>2/ check if all the GUI and make sure it works</li> <li>3/ check if the application closes without errors</li> </ul>

<b>Open Issues:</b>	
---------------------	--

<b>Use Case Number:</b>	2
<b>Use Case Name:</b>	<i>Take Photo</i>
<b>Participating Actors:</b>	main User (patient)
<b>Goal:</b>	Allow the User to take a photo, and then have the application store the photo with a storage method.
<b>Trigger:</b>	The User chooses <u>take photo</u> option.
<b>Precondition:</b>	User knows the target area; there is a camera.
<b>Postcondition:</b>	On success, all relevant information of the captured photo is stored.
<b>Basic Flow</b>	1/ User takes a photo 2/ Use case 3: Timestamp A Photo 3/ System displays the taken photo 4/ System prompts the User to enter the name of a new <i>skin condition</i> or to choose some existing <i>skin condition</i> that the taken photo refers to 5/ User enters the name of a new <i>skin condition</i> or choose some existing <i>skin condition</i> 6/ System prompts the User to enter the name of a new <i>group</i> or to choose some existing <i>group</i> for the taken photo 7/ User enters the name of a new <i>group</i> or choose some existing <i>group</i> 8/ System prompts the User to enter Annotation to be associated with the taken photo 9/ User enters some annotation about the photo 10/ System saves the taken photo
<b>Exceptions</b>	1/ If the camera malfunctions 1.1/ System displays an error 1.2/ System returns to the previous activity 8/ If the System can't save the taken photo due to perhaps insufficient amount of storage 8.1/ System displays an error 8.2/ System prompts the User to clear some memory or to discard the

	<p>taken photo</p> <p>8.3/ If the User chooses <u>clear some memory</u> return to step 3</p> <p>8.4/ If the User chooses <u>discard the taken photo</u> return to the previous activity</p>
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	Timestamp A Photo, Tag a photo, add skin conditions to photo
<b>Extends:</b>	
<b>Related Artifacts:</b>	<p>As a user, I want pictures I take to be stored so that I may recall them and view them later.</p> <p>Camera</p>
<b>Notes:</b>	<p><b><u>-Tests for use case 2-</u></b></p> <p>1/ verify the application can capture an image with its camera</p> <p>2/ verify the photo can be displayed on screen</p> <p>3/ test to see if all user interaction with the User Interface works correctly such as user input, and the user selections on the skin condition, and group options.</p> <p>4/ verify that the picture can be stored</p> <p>5/ verify that the picture can be reloaded after being stored</p> <p>6/ verify how the program handles running out of memory to store a new photo</p>
<b>Open Issues:</b>	Can we implement a better looking camera?

<b>Use Case Number:</b>	3
<b>Use Case Name:</b>	<i>Timestamp A Photo</i>
<b>Participating Actors:</b>	This is done automatically by the system

<b>Goal:</b>	Timestamp A Photo
<b>Trigger:</b>	The User has taken a photo
<b>Precondition:</b>	The photo is taken successfully by the camera
<b>Postcondition:</b>	On success, System automatically records the time stamp of the taken photo
<b>Basic Flow</b>	1/ Use case 2 2/ System automatically records the time stamp on and add to the taken photo
<b>Exceptions</b>	2/ If the System cannot add a time stamp 2.1/ System notifythe user that it cannot record the time stamp. 2.2/ System does not record time stamp of the photo.
<b>Qualities:</b>	Time stamp is recoded with the following format DD/MM/YYYY HH:MM
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want all pictures I take to be timestamped so that I can measure the time between photos and know when an event occurred.
<b>Notes:</b>	<p>This use case is an activity inside of the “take a photo” class, it will automatically add a timestamp to any photos taken by the user that they wish to save to memory, there are many simple java methods to implement this task.</p> <p><b><u>-Tests for use case 3-</u></b></p> <p>1/ check that stored photos have a time stamp 2/ check that they have the correct time stamp</p>

<b>Open Issues:</b>	
---------------------	--

<b>Use Case Number:</b>	4
<b>Use Case Name:</b>	<i>View Photo</i>
<b>Participating Actors:</b>	main User (patient)
<b>Goal:</b>	Allow the User to review a taken photo
<b>Trigger:</b>	The User chooses <u>view photo</u> option
<b>Precondition:</b>	User knows which photo he/she wishes to review
<b>Postcondition:</b>	On success, all relevant information of the captured photo ( <i>skin condition, body part and timestamp</i> ) is displayed
<b>Basic Flow</b>	1/ System prompts the User to choose a photo to view 2/ User chooses a photo to view 3/ System displays the chosen photo 4/ System gives option to view photo annotation and other misc information about the photo
<b>Exceptions</b>	3/ If the System cannot display the chosen photo because the photo may not exist or it may be broken 3.1/ System notify the user that it cannot display the chosen image 3.2/ System returns to step 1
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	

<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want to be able to review any photos I have taken and view them so that I can see any changes.
<b>Notes:</b>	<u><b>-Tests for use case 4-</b></u> 1/ verify the user may view a photo
<b>Open Issues:</b>	

<b>Use Case Number:</b>	5
<b>Use Case Name:</b>	<i><b>Compare Two Photos</b></i>
<b>Participating Actors:</b>	main User
<b>Goal:</b>	View two photos simultaneously to observe any change of some skin condition
<b>Trigger:</b>	The User chooses <u><i>compare multiple photos</i></u> option
<b>Precondition:</b>	The User knows up front which photos he/she wishes to compare
<b>Postcondition:</b>	On success, System displays two photos simultaneously
<b>Basic Flow</b>	1/ User choose to view all photos OR to choose a <i>skin condition</i> in a list to show all photos of the chosen <i>skin condition</i> ( <i>use case View skin condition</i> ) OR to choose a <i>group</i> in a list to show all photos of the chosen <i>group</i> ( <i>use case View group</i> ) 2/ System prompts User to choose two photos for comparison 3/ User choose to compare photos 4/ System shows two photos simultaneously by default setting 5/ User to choose type of view ( side-by-side or overlay) 6/ System shows photos side-by-side or overlay depending on User's choice

<b>Exceptions</b>	2/ User chooses one or more than two photos 2.1/ Ask the user to choose two photos to compare
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	Only 2 photos can be viewed simultaneously
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want a method to organize these pictures so that I can view pictures of the same bit of skin over time and see if there is any progression or growth. For instance I might want to organize some photos that are related to my "gross mole that occurs on my right hand".
<b>Notes:</b>	<b><u>-Tests for use case 5-</u></b> 1/ verify the application can display two photos properly 2/ verify the program acts appropriately when the user only has 0 images, and when the user only has 1 image 3/ verify the user can choose two photos 4/ verify the view works properly for displaying the two photos
<b>Open Issues:</b>	

<b>Use Case Number:</b>	6
<b>Use Case Name:</b>	<i>Tag a photo</i>
<b>Participating Actors:</b>	main User (patient)
<b>Goal:</b>	Add a photo to an existing group or a new tag
<b>Trigger:</b>	User selects <u>tag a photo</u>



<b>Precondition:</b>	User knows which photo to add to some tag
<b>Postcondition:</b>	On success, the chosen photo is added to some tag
<b>Basic Flow</b>	1/ System prompts the User to choose a photo 2/ User selects a photo 3/ System prompts User to assign some existing <i>tag</i> or create a new <i>tag</i> for the chosen photo 4/ The user chooses existing <i>tag</i> or the newly created <i>tag</i> for the taken photo
<b>Exceptions</b>	3/ If there is an error when making a new tag 3.1/ System shows the error message 3.2/ Return to the previous activity
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user if I have an organizational tag or group, I want to be able to add photos to that group so that I can further track the progression of something relevant to that group. # This use case happens when the user hasn't placed some photo in a <i>group</i> when the photo was taken
<b>Notes:</b>	<u><b>-Tests for use case 6-</b></u> 1/ verify that a photo can be added to a tag 2/ handle the null case, where the user tries to add null to a tag?
<b>Open Issues:</b>	

<b>Use Case Number:</b>	7
-------------------------	---

<b>Use Case Name:</b>	<i>Create Transparent Layer</i>
<b>Participating Actors:</b>	User
<b>Goal:</b>	Use some photo as a transparent layer to help user take a new photo in the same position as the layer
<b>Trigger:</b>	User selects <i>transparent layer</i> option inside the <i>taking photo</i> view
<b>Precondition:</b>	User knows the photo to make it a transparent layer
<b>Postcondition:</b>	On success, a <i>transparent layer</i> appears on a screen
<b>Basic Flow</b>	1/ System prompts the User to choose a photo 2/ User selects a photo 3/ System displays the <i>transparent layer</i> on the screen
<b>Exceptions</b>	3/ If there is an error when making a new group 3.1/ System shows the error message 3.2/ Return to the previous activity
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	Take Photo
<b>Related Artifacts:</b>	As a user, I want some method of helping me take consistent photos, so that when I show the doctor any progression such as the growth of a mole is evident.

<b>Notes:</b>	<p>This use case happens when the user wishes to take a new photo based on some existing photo of the same skin condition.</p> <p><u><b>-Tests for use case 7-</b></u></p> <p>1/ verify that the layer works</p> <p>2/ verify a photo can still be capture with the layer</p> <p>3/ handle the case when no photos are present for the user to use as a layer</p>
<b>Open Issues:</b>	This use case requires knowledge of image processing, thus it will be added later if possible.

<b>Use Case Number:</b>	8
<b>Use Case Name:</b>	<i>Create Reminder</i>
<b>Participating Actors:</b>	User
<b>Goal:</b>	Remind user of taking photo of some skin condition regularly
<b>Trigger:</b>	User choose <i>set reminder</i> option
<b>Precondition:</b>	User knows the skin condition and the group
<b>Postcondition:</b>	The System sets the alarm to remind User of taking photo and specify the skin condition and group automatically based on the info of the alarm
<b>Basic Flow</b>	<p>1/ System prompts the User to choose the skin condition in a list</p> <p>2/ User chooses some skin condition</p> <p>3/ System prompts the User to choose the group in a list</p> <p>4/ User chooses some group</p> <p>5/ System prompts the User to enter other information of the alarm: time, repetition, name, additional info.</p> <p>6/ User change the default value of these information if needed</p> <p>7/ System sets the reminder</p>
<b>Exceptions</b>	<p>7/ If System cannot set the reminder</p> <p>7.1/ System displays an error message</p> <p>7.2/ Return to previous activity</p>
<b>Qualities:</b>	
<b>Functionality:</b>	

<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want some method of helping me take consistent photos, so that when I show the doctor any progression such as the growth of a mole is evident.
<b>Notes:</b>	<u><b>-Tests for use case 8-</b></u> 1/ verify a reminder can be created
<b>Open Issues:</b>	

<b>Use Case Number:</b>	9
<b>Use Case Name:</b>	<i><b>Retake Photo</b></i>
<b>Participating Actors:</b>	User
<b>Goal:</b>	Retake an erroneous photo
<b>Trigger:</b>	User choose to retake a photo
<b>Precondition:</b>	The user has and able to modify the photo he or she want to retake
<b>Postcondition:</b>	The old unwanted photo is replaced by a new photo
<b>Basic Flow</b>	1/ User choose the option to retake that photo 2/ System invokes the function in use case 2 3/ New photo overwrites the old photo
<b>Exceptions</b>	3/ If the new photo cannot overwrite the old photo 3.1/ System displays an error message 3.2/ The user will then have to go and delete the previously taken photo

<b>Qualities:</b>	The process of retaking a photo is simple and quick.
<b>Functionality:</b>	
<b>Constraints:</b>	The new photo must reuse any appropriate information attached in the old photo.
<b>Includes:</b>	Take Photo
<b>Extends:</b>	
<b>Related Artifacts:</b>	# As a user, I should be able to retake photos I am taking, if I fail to take a photo # I want to correct it, so that I do not have erroneous photos. --Camera
<b>Notes:</b>	<u><b>-Tests for use case 9-</b></u> 1/ verify the user is able to use the retake photo option 2/ verify that the previous photo has been removed from the application 3/ verify the new photo is able to be stored by the application
<b>Open Issues:</b>	

<b>Use Case Number:</b>	10
<b>Use Case Name:</b>	<i>View Skin Conditions</i>
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of skin conditions
<b>Trigger:</b>	The User chooses <u>view list of skin conditions</u>
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of skin conditions or an empty list if there is no skin condition
<b>Basic Flow</b>	1/ System displays the list of skin conditions
<b>Exceptions</b>	1/ If System cannot display the list of skin conditions 1.1/ System displays an error message 1.2/ Return to previous activity

<b>Qualities:</b>	The list must be easy to read.
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	View Multiple Photos
<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	<p><b><u>-Tests for use case 10-</u></b></p> <p>1/ verify all the skin conditions are displayed correctly</p> <p>2/ verify the case where no skin conditions exist is handled correctly</p>
<b>Open Issues:</b>	

<b>Use Case Number:</b>	11
<b>Use Case Name:</b>	<i>View Groups</i>
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of groups
<b>Trigger:</b>	The User chooses <u>view list of tags</u>
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of groups or an empty list if there is no tag
<b>Basic Flow</b>	1/ System displays the list of tags
<b>Exceptions</b>	1/ If System cannot display the list of tags 1.1/ System displays an error message 1.2/ Return to previous activity



<b>Qualities:</b>	The list must be easy to read
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	View Multiple Photos
<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	<p><b><u>-Tests for use case 11-</u></b></p> <p>1/ verify all the appropriate tags are displayed</p> <p>2/ verify the case is handled where there are no tags</p>
<b>Open Issues:</b>	

<b>Use Case Number:</b>	12
<b>Use Case Name:</b>	<i>View Photos in a Skin Condition</i>
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of photos in a chosen skin condition
<b>Trigger:</b>	The User clicks on a skin condition in the list of skin conditions
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of photos in the chosen skin condition
<b>Basic Flow</b>	1/ User chooses a skin condition in the list 2/ System displays a list of photos in the chosen skin condition 3/ User click on a photo to view it bigger
<b>Exceptions</b>	The list of photos must be loaded fast and correctly

<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	View Skin Conditions
<b>Related Artifacts:</b>	
<b>Notes:</b>	<u><b>Tests for use case 12-</b></u> 1/ verify the case is handled when no skin condition is selected 2/ check to see if all the photos with that skin condition are displayed
<b>Open Issues:</b>	

<b>Use Case Number:</b>	13
<b>Use Case Name:</b>	<i>View Photos in Tag</i>
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of photos in a chosen tag
<b>Trigger:</b>	The User clicks on a group in the list of tags
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of photos in the chosen tag
<b>Basic Flow</b>	1/ User chooses a group in the list of tags 2/ System displays a list of photos in the chosen tag 3/ User click on a photo to view it bigger
<b>Exceptions</b>	

<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	View tags
<b>Related Artifacts:</b>	
<b>Notes:</b>	<p><b><u>-Tests for use case 13-</u></b></p> <p>1/ verify the case where no tags are selected</p> <p>2/ verify that the correct photos are displayed corresponding to the tag</p>
<b>Open Issues:</b>	

<b>Use Case Number:</b>	14
<b>Use Case Name:</b>	<i>Edit Photo Info</i>
<b>Participating Actors:</b>	main user (patient)
<b>Goal:</b>	Edit information of a photo
<b>Trigger:</b>	User chooses the <u>edit</u> option in the menu
<b>Precondition:</b>	User know what information he/she want to edit
<b>Postcondition:</b>	On success, all relevant information of the chosen photo is updated
<b>Basic Flow</b>	1/ System has displayed info of the chosen photo at the end of use case <u>View photo</u> 2/ User chooses some piece of information of the entry to edit (time stamp is not editable) (optional) 3/ User confirms new information 4/ System saves and updates the entry

<b>Exceptions</b>	3/ If the new information is not valid 3.1/ System displays an error 3.2/ System returns to step 2
<b>Qualities:</b>	System assists user in editing by providing appropriate keys
<b>Functionality:</b>	
<b>Constraints:</b>	Use DD/MM/YYYY format for the date of a log entry; a short mnemonic name is no more than 100 characters
<b>Includes:</b>	
<b>Extends:</b>	View Photo
<b>Related Artifacts:</b>	
<b>Notes:</b>	<u><b>-Tests for use case 14-</b></u> 1/ verify the photo is saved properly after being edited 2/ verify the new fields are updated correctly 3/ handle the case where the new fields are null

<b>Open Issues:</b>	
---------------------	--

<b>Use Case Number:</b>	15
<b>Use Case Name:</b>	<i>Delete Photo</i>
<b>Participating Actors:</b>	main user (patient)
<b>Goal:</b>	Delete an existing photo
<b>Trigger:</b>	User chooses the <i>Delete</i> option in the menu after choosing some photo
<b>Precondition:</b>	User knows what photo he/she want to delete
<b>Postcondition:</b>	On success, delete the photo



<b>Basic Flow</b>	1/ System has displayed the chosen photo at the end of <i>View photo</i> use case 2/ User chooses “delete” option to delete the photo 3/ System prompts the user to confirm action 4/ Confirmation of delete action
<b>Exceptions</b>	4/ if System cannot delete the photo 4.1/ System shows an error message 4.2/ Return to step 1
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	

<b>Notes:</b>	<u><b>-Tests for use case 15-</b></u> 1/ verify that the photo is removed, and all other data relating to that photo is also removed
<b>Open Issues:</b>	

<b>Use Case Number:</b>	16
<b>Use Case Name:</b>	Delete Tag/Skin Condition
<b>Participating Actors:</b>	main user (patient)
<b>Goal:</b>	Delete some tag or some skin condition
<b>Trigger:</b>	User chooses the <u>Delete</u> option in the menu after the System shows a list of tags or a list of skin conditions
<b>Precondition:</b>	

<b>Postcondition:</b>	On success, remove the chosen tag or skin condition out of the information of all photos that belong to the tag or skin condition
<b>Basic Flow</b>	1/ System has displayed the list of groups or the list of skin conditions at the end of <u>View tags</u> or <u>View skin conditions</u> 2/ User chooses <u>delete</u> option in the menu 3/ System prompts User to choose tags or skin conditions to delete 4/ User chooses tags or skin conditions to delete 5/ System prompts User to confirm deletion 6/ User confirms deletion 7/ System identifies all photos that belong to the deleted tag or deleted skin condition and remove the group or skin condition out of the information of those photos; remove the tag or skin condition
<b>Exceptions</b>	7/ if System cannot modify info of the photos 7.1/ System shows an error message 7.2/ Not delete the tag or skin condition 7.3/ Return to previous activity
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	

<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	<p><u><b>-Tests for use case 16-</b></u></p> <p>1/ verify that the tag/skin condition is deleted</p> <p>2/ verify that all photos who were in the tag or had the skin condition tags now no longer have the tags</p>
<b>Open Issues:</b>	

<b>Use Case Number:</b>	17
<b>Use Case Name:</b>	View Multiple Photos
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of photos

<b>Trigger:</b>	The User chooses a criteria to view a list of photos
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of photos or an empty list if there is no photo meet the criteria
<b>Basic Flow</b>	1/ System displays a list of photos meet the criteria 2/ User click on a photo to view it bigger
<b>Exceptions</b>	
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	

<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	<u><i>-Tests for use case 17-</i></u> 1/ verify that the list of photos is displayed correctly 2/ verify that the list of photos for 0 photos doesn't create errors
<b>Open Issues:</b>	

<b>Use Case Number:</b>	18
<b>Use Case Name:</b>	<i>Password encryption</i>
<b>Participating Actors:</b>	Main user
<b>Goal:</b>	Allow returned user to log in with his/her password

<b>Trigger:</b>	This user clicks the app from the android app menu list
<b>Precondition:</b>	None
<b>Postcondition:</b>	Allow user to log in if the password is correct, reject otherwise
<b>Basic Flow</b>	1/ User click application icon 2/ System prompts user for the password 3/ User enters password 4/ System checks the password 5/ System starts the application if password is correct. Otherwise, display a wrong password message. Return to step 2
<b>Exceptions</b>	2.1/ application fails to launch Display error message
<b>Qualities:</b>	System responds within 1s
<b>Functionality:</b>	
<b>Constraints:</b>	

<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I am worried that if my phone gets stolen people will see my photos. I want my photos secured via a password
<b>Notes:</b>	Implement encryption of the folder containing photos
<b>Object Oriented Analysis</b>	

<b>Use Case Number:</b>	19
<b>Use Case Name:</b>	<i>Reminder</i>
<b>Participating Actors:</b>	Main user



<b>Goal:</b>	Remind the user of taking consistent photos to keep track of his/her skin conditions.
<b>Trigger:</b>	This user chooses “reminder” option on the main screen of the application
<b>Precondition:</b>	None
<b>Postcondition:</b>	Allow user to create alarms that go off at set-up times.
<b>Basic Flow</b>	1/ User chooses “reminder” option 2/ System shows a list of alarms created 3/ User chooses to create a new alarm or edit some existing one 4/ System creates a new alarm or modifies the existing one
<b>Exceptions</b>	2/ Fails to load the list of alarms 2.1/ Display an error message 2.2/ Return to the main screen
<b>Qualities:</b>	
<b>Functionality:</b>	The alarm goes of at a specific point of time that is set by the user. If the time set up already passes the current time, the alarm will go off immediately. One-time or repeating alarm types are supported.
<b>Constraints:</b>	

<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want notifications about when I should take photos again. I'm forgetful and being told every so often to take new photos would be helpful. Certain photos I don't have to retake again.
<b>Notes:</b>	The alarm is not bound to any particular photos. The user needs to annotate the alarm to know which skin condition he/she would like to observe by taking consistent photos. The alarm name is stored in a database
<b>Open Issues:</b>	

<b>Use Case Number:</b>	20
<b>Use Case Name:</b>	<i>Emailing photos</i>
<b>Participating Actors:</b>	Main user

<b>Goal:</b>	Email a set of photos to the patient's physician
<b>Trigger:</b>	This user chooses "email photos" option to email a set of chosen photos to his/her physician
<b>Precondition:</b>	There is a list of photos being displayed on the screen
<b>Postcondition:</b>	Successfully email the physician a set of chosen photos
<b>Basic Flow</b>	1/ User chooses "email photos" option 2/ System prompts user to choose multiple photos within the list 3/ User choose a list of photos 4/ User confirms his/her chosen photo set 5/ System displays pop-up window to ask for the user's confirmation of sending photos via email 6/ Upon confirmation of the user, system prompts the user to enter the physician's email 7/ System sends the chosen set of photos to the provided email 8/ System notifies the user of successful or failed email
<b>Exceptions</b>	
<b>Qualities:</b>	

<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want to email a set of photos to my physician
<b>Notes:</b>	Checking the network connection before sending photos via email to provide descriptive feedback
<b>Open Issues:</b>	

<b>Use Case Number:</b>	21
-------------------------	----

<b>Use Case Name:</b>	<i>Annotate images</i>
<b>Participating Actors:</b>	Main user
<b>Goal:</b>	Annotate images to provide more relevant information about the images
<b>Trigger:</b>	This user chooses “add annotation” option when the newly created photo is displayed on the screen right after it is taken
<b>Precondition:</b>	User takes a photo
<b>Postcondition:</b>	Successfully annotate the photo
<b>Basic Flow</b>	1/ User chooses “add annotation” option 2/ System prompts user to input some information 3/ User inputs an annotation and confirms his/her input 4/ System attaches the annotation to the photo
<b>Exceptions</b>	
<b>Qualities:</b>	

<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user I would like to annotate images with annotations to describe the scenario that I made this photo
<b>Notes:</b>	The annotation is stored in the database. It will be loaded when the user chooses “show annotation” option when viewing the photo
<b>Open Issues:</b>	

<b>Use Case Number:</b>	22
-------------------------	----

<b>Use Case Name:</b>	<i>Tag photos</i>
<b>Participating Actors:</b>	Main user
<b>Goal:</b>	Tag the images to provide queries images by tag
<b>Trigger:</b>	The user chooses some existing tag to tag the photo to or create a new tag for the photo
<b>Precondition:</b>	User has just taken a photo or some tag is chosen
<b>Postcondition:</b>	Successfully tag the photo
<b>Basic Flow</b>	<p>When a photo is just taken</p> <ol style="list-style-type: none"> <li>1/ System loads the list of tags</li> <li>2/ User chooses some existing tags to tag the photo</li> <li>3/ User confirms his/her chosen tags</li> <li>4/ System tags photo with chosen tags</li> </ol> <p>The target tag can be chosen from the list of existing tags</p> <ol style="list-style-type: none"> <li>1/ User chooses “connect photos” option to tag a chosen photo to the current selected tag</li> <li>2/ System loads a list of photos that are currently not attached to the tag</li> <li>3/ User chooses some photo(s) to tag</li> <li>4/ System tags the chosen photos</li> </ol>

<b>Exceptions</b>	
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user I would like to tag the images so I can query the image by tag
<b>Notes:</b>	
<b>Open Issues:</b>	



