

<b>Use Case Number:</b>	1
<b>Use Case Name:</b>	Start Skin Viewer application
<b>Participating Actors:</b>	main user (patient)
<b>Goal:</b>	initialize storage methods/check for previously stored data, open up the application for the user
<b>Trigger:</b>	this user clicks the app from the android app menu list
<b>Precondition:</b>	none
<b>Postcondition:</b>	The application is now at start up
<b>Basic Flow</b>	1/ User click application icon 2/ Application is launched
<b>Exceptions</b>	2.1/ application fails to launch display error message
<b>Qualities:</b>	System responds within 1s
<b>Functionality:</b>	<p>The application has:</p> <ul style="list-style-type: none"> <li>search/add/delete/view/edit</li> <li>---&gt; 1/ search: photo, skin condition, group</li> <li>---&gt; 2/ create: a photo, a skin condition(entry)</li> <li>---&gt; 3/ add: entries to group, photos to entries</li> <li>---&gt; 4/ delete: photos, entries, groups</li> <li>---&gt; 5/ view: photos, skin conditions, groups,</li> <li>---&gt; 6/ edit: photos, skin conditions, groups</li> </ul>
<b>Constraints:</b>	The UI looks aesthetically pleasing
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want to keep a visual log of any skin conditions I might . The purpose could be to gather a personal history of a medical condition as to provide evidence to a physician.
<b>Notes:</b>	We'll add details into each function later

<b>Open Issues:</b>	
---------------------	--

<b>Use Case Number:</b>	2
<b>Use Case Name:</b>	Take Photo
<b>Participating Actors:</b>	main User (patient)
<b>Goal:</b>	Allow the User to take a photo, and then have the application store the photo with a storage method.
<b>Trigger:</b>	The User chooses <u>take photo</u> option.
<b>Precondition:</b>	User knows the target area; there is a camera.
<b>Postcondition:</b>	On success, all relevant information of the captured photo is stored.
<b>Basic Flow</b>	1/ User takes a photo 2/ Use case 3: Timestamp A Photo 3/ System displays the taken photo 4/ System prompts the User to enter the name of a new <i>skin condition</i> or to choose some existing <i>skin condition</i> that the taken photo refers to 5/ User enters the name of a new <i>skin condition</i> or choose some existing <i>skin condition</i> 6/ System prompts the User to enter the name of a new <i>group</i> or to choose some existing <i>group</i> for the taken photo 7/ User enters the name of a new <i>group</i> or choose some existing <i>group</i> 8/ System saves the taken photo
<b>Exceptions</b>	1/ If the camera malfunctions 1.1/ System displays an error 1.2/ System returns to the previous activity 8/ If the System can't save the taken photo due to perhaps insufficient amount of storage 8.1/ System displays an error 8.2/ System prompts the User to clear some memory or to discard the taken photo 8.3/ If the User chooses <u>clear some memory</u> return to step 3 8.4/ If the User chooses <u>discard the taken photo</u> return to the previous activity
<b>Qualities:</b>	The system must notify the user the image has been stored successfully.
<b>Functionality:</b>	
<b>Constraints:</b>	

<b>Includes:</b>	Timestamp A Photo
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want pictures I take to be stored so that I may recall them and view them later. Camera
<b>Notes:</b>	
<b>Open Issues:</b>	How should a photo be named?

<b>Use Case Number:</b>	3
<b>Use Case Name:</b>	Timestamp A Photo
<b>Participating Actors:</b>	
<b>Goal:</b>	Timestamp A Photo
<b>Trigger:</b>	The User has taken a photo
<b>Precondition:</b>	The photo is taken successfully by the camera
<b>Postcondition:</b>	On success, System automatically records the time stamp of the taken photo
<b>Basic Flow</b>	1/ Use case 2 2/ System automatically records the time stamp on and add to the taken photo
<b>Exceptions</b>	2/ If the System cannot add a time stamp 2.1/ System notifythe user that it cannot record the time stamp. 2.2/ System does not record time stamp of the photo.
<b>Qualities:</b>	Time stamp is recoded with the following format HH:MM DD/MM/YYYY
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want all pictures I take to be timestamped so that I can measure the time between photos and know when an event occurred.
<b>Notes:</b>	This use case is an activity inside of the “take a photo” class, it will automatically add a timestamp to any photos taken by the user that they wish to save to

	memory, there are many simple java methods to implement this task.
<b>Open Issues:</b>	

<b>Use Case Number:</b>	4
<b>Use Case Name:</b>	View Photo
<b>Participating Actors:</b>	main User (patient)
<b>Goal:</b>	Allow the User to review a taken photo
<b>Trigger:</b>	The User chooses <u>view photo</u> option
<b>Precondition:</b>	User knows which photo he/she wishes to review
<b>Postcondition:</b>	On success, all relevant information of the captured photo ( <i>skin condition, body part and timestamp</i> ) is displayed
<b>Basic Flow</b>	1/ System prompts the User to choose a photo to view 2/ User chooses a photo to view 3/ System displays the chosen photo
<b>Exceptions</b>	3/ If the System cannot display the chosen photo because the photo may not exist or it may be broken 3.1/ System notify the user that it cannot display the chosen image 3.2/ System returns to step 1
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want to be able to review any photos I have taken and view them so that I can see any changes.

<b>Notes:</b>	
<b>Open Issues:</b>	



<b>Use Case Number:</b>	5
<b>Use Case Name:</b>	Compare Two Photos
<b>Participating Actors:</b>	main User
<b>Goal:</b>	View two photos simultaneously to observe any change of some skin condition
<b>Trigger:</b>	The User chooses <i>compare multiple photos</i> option
<b>Precondition:</b>	The User knows up front which photos he/she wishes to compare
<b>Postcondition:</b>	On success, System displays two photos simultaneously
<b>Basic Flow</b>	1/ User choose to view all photos OR to choose a <i>skin condition</i> in a list to show all photos of the chosen <i>skin condition</i> ( <i>use case View skin condition</i> ) OR to choose a <i>group</i> in a list to show all photos of the chosen <i>group</i> ( <i>use case View group</i> ) 2/ System prompts User to choose two photos for comparison 3/ User choose to compare photos 4/ System shows two photos simultaneously by default setting 5/ User to choose type of view ( side-by-side or overlay) 6/ System shows photos side-by-side or overlay depending on User's choice
<b>Exceptions</b>	2/ User chooses one or more than two photos 2.1/ Ask the user to choose two photos to compare
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	Only 2 photos can be viewed simultaneously
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related</b>	As a user, I want a method to organize these pictures so that I can view pictures

<b>Artifacts:</b>	of the same bit of skin over time and see if there is any progression or growth. For instance I might want to organize some photos that are related to my "gross mole that occurs on my right hand".
<b>Notes:</b>	We need some method to list the photos, and we need some method to order the list, like a SQL ORDER BY command.
<b>Open Issues:</b>	Comparison of more than 2 photos will be considered.

<b>Use Case Number:</b>	6
<b>Use Case Name:</b>	Add Photo to a Group
<b>Participating Actors:</b>	main User (patient)
<b>Goal:</b>	Add a photo to an existing group or a new group
<b>Trigger:</b>	User selects <u>add photo to group</u>
<b>Precondition:</b>	User knows which photo to add to some group
<b>Postcondition:</b>	On success, the chosen photo is added to some group
<b>Basic Flow</b>	1/ System prompts the User to choose a photo 2/ User selects a photo 3/ System prompts User to assign some existing <i>group</i> or create a new <i>group</i> for the chosen photo 4/ The user chooses existing <i>group</i> or the newly created <i>group</i> for the taken photo
<b>Exceptions</b>	3/ If there is an error when making a new group 3.1/ System shows the error message 3.2/ Return to the previous activity
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user if I have an organizational tag or group, I want to be able to add photos to that group so that I can further track the progression of something relevant to that group. # This use case happens when the user hasn't placed some photo in a <i>group</i> when the photo was taken

<b>Notes:</b>	
<b>Open Issues:</b>	

<b>Use Case Number:</b>	7
<b>Use Case Name:</b>	Create Transparent Layer
<b>Participating Actors:</b>	User
<b>Goal:</b>	Use some photo as a transparent layer to help user take a new photo in the same position as the layer
<b>Trigger:</b>	User selects <i>transparent layer</i> option inside the <i>taking photo</i> view
<b>Precondition:</b>	User knows the photo to make it a transparent layer
<b>Postcondition:</b>	On success, a <i>transparent layer</i> appears on a screen
<b>Basic Flow</b>	1/ System prompts the User to choose a photo 2/ User selects a photo 3/ System displays the <i>transparent layer</i> on the screen
<b>Exceptions</b>	3/ If there is an error when making a new group 3.1/ System shows the error message 3.2/ Return to the previous activity
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	Take Photo
<b>Related Artifacts:</b>	As a user, I want some method of helping me take consistent photos, so that when I show the doctor any progression such as the growth of a mole is evident.
<b>Notes:</b>	This use case happens when the user wishes to take a new photo based on some existing photo of the same skin condition

<b>Open Issues:</b>	This use case requires knowledge of image processing, thus it will be added later if possible.
---------------------	--

<b>Use Case Number:</b>	8
<b>Use Case Name:</b>	Create Reminder
<b>Participating Actors:</b>	User
<b>Goal:</b>	Remind user of taking photo of some skin condition regularly
<b>Trigger:</b>	User choose <i>set reminder</i> option
<b>Precondition:</b>	User knows the skin condition and the group
<b>Postcondition:</b>	The System sets the alarm to remind User of taking photo and specify the skin condition and group automatically based on the info of the alarm
<b>Basic Flow</b>	1/ System prompts the User to choose the skin condition in a list 2/ User chooses some skin condition 3/ System prompts the User to choose the group in a list 4/ User chooses some group 5/ System prompts the User to enter other information of the alarm: time, repetition, name, additional info. 6/ User change the default value of these information if needed 7/ System sets the reminder
<b>Exceptions</b>	7/ If System cannot set the reminder 7.1/ System displays an error message 7.2/ Return to previous activity
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	As a user, I want some method of helping me take consistent photos, so that when I show the doctor any progression such as the growth of a mole is evident.
<b>Notes:</b>	
<b>Open Issues:</b>	

<b>Use Case Number:</b>	9
<b>Use Case Name:</b>	Retake Photo
<b>Participating Actors:</b>	User
<b>Goal:</b>	Retake an erroneous photo
<b>Trigger:</b>	User choose to retake a photo
<b>Precondition:</b>	The user has and able to modify the photo he or she want to retake
<b>Postcondition:</b>	The old unwanted photo is replaced by a new photo
<b>Basic Flow</b>	1/ User choose the option to retake that photo 2/ System invokes the function in use case 2 3/ New photo overwrites the old photo
<b>Exceptions</b>	3/ If the new photo cannot overwrite the old photo 3.1/ System displays an error message 3.2/ The user will then have to go and delete the previously taken photo
<b>Qualities:</b>	The process of retaking a photo is simple and quick.
<b>Functionality:</b>	
<b>Constraints:</b>	The new photo must reuses any appropriate information attached in the old photo.
<b>Includes:</b>	Take Photo
<b>Extends:</b>	
<b>Related Artifacts:</b>	# As a user, I should be able to retake photos I am taking, if I fail to take a photo # I want to correct it, so that I do not have erroneous photos. --Camera
<b>Notes:</b>	
<b>Open Issues:</b>	



<b>Use Case Number:</b>	10
<b>Use Case Name:</b>	View Skin Conditions
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of skin conditions
<b>Trigger:</b>	The User chooses <i><u>view list of skin conditions</u></i>
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of skin conditions or an empty list if there is no skin condition
<b>Basic Flow</b>	1/ System displays the list of skin conditions
<b>Exceptions</b>	1/ If System cannot display the list of skin conditions 1.1/ System displays an error message 1.2/ Return to previous activity
<b>Qualities:</b>	The list must be easy to read.
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	View Multiple Photos
<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	
<b>Open Issues:</b>	

<b>Use Case Number:</b>	11
<b>Use Case Name:</b>	View Groups
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of groups
<b>Trigger:</b>	The User chooses <i><u>view list of groups</u></i>
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of groups or an empty list if there is no group
<b>Basic Flow</b>	1/ System displays the list of groups
<b>Exceptions</b>	1/ If System cannot display the list of groups 1.1/ System displays an error message 1.2/ Return to previous activity
<b>Qualities:</b>	The list must be easy to read
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	View Multiple Photos
<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	
<b>Open Issues:</b>	

<b>Use Case Number:</b>	12
<b>Use Case Name:</b>	View Photos in a Skin Condition
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of photos in a chosen skin condition
<b>Trigger:</b>	The User clicks on a skin condition in the list of skin conditions
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of photos in the chosen skin condition
<b>Basic Flow</b>	1/ User chooses a skin condition in the list 2/ System displays a list of photos in the chosen skin condition 3/ User click on a photo to view it bigger
<b>Exceptions</b>	The list of photos must be loaded fast and correctly
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	View Skin Conditions
<b>Related Artifacts:</b>	
<b>Notes:</b>	
<b>Open Issues:</b>	

<b>Use Case Number:</b>	13
<b>Use Case Name:</b>	View Photos in Group
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of photos in a chosen group
<b>Trigger:</b>	The User clicks on a group in the list of groups
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of photos in the chosen group
<b>Basic Flow</b>	1/ User chooses a group in the list of groups 2/ System displays a list of photos in the chosen group 3/ User click on a photo to view it bigger
<b>Exceptions</b>	
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	View Groups
<b>Related Artifacts:</b>	
<b>Notes:</b>	
<b>Open Issues:</b>	

<b>Use Case Number:</b>	14
<b>Use Case Name:</b>	Edit Photo Info
<b>Participating Actors:</b>	main user (patient)
<b>Goal:</b>	Edit information of a photo
<b>Trigger:</b>	User chooses the <i><u>edit</u></i> option in the menu
<b>Precondition:</b>	User know what information he/she want to edit
<b>Postcondition:</b>	On success, all relevant information of the chosen photo is updated
<b>Basic Flow</b>	1/ System has displayed info of the chosen photo at the end of use case <i><u>View photo</u></i> 2/ User chooses some piece of information of the entry to edit (time stamp is not editable) (optional) 3/ User confirms new information 4/ System saves and updates the entry
<b>Exceptions</b>	3/ If the new information is not valid 3.1/ System displays an error 3.2/ System returns to step 2
<b>Qualities:</b>	System assists user in editing by providing appropriate keys
<b>Functionality:</b>	
<b>Constraints:</b>	Use DD/MM/YYYY format for the date of a log entry; a short mnemonic name is no more than 100 characters
<b>Includes:</b>	
<b>Extends:</b>	View Photo
<b>Related Artifacts:</b>	
<b>Notes:</b>	
<b>Open Issues:</b>	

<b>Use Case Number:</b>	15
<b>Use Case Name:</b>	Delete Photo
<b>Participating Actors:</b>	main user (patient)
<b>Goal:</b>	Delete an existing photo
<b>Trigger:</b>	User chooses the <i>Delete</i> option in the menu after choosing some photo
<b>Precondition:</b>	User knows what photo he/she want to delete
<b>Postcondition:</b>	On success, delete the photo
<b>Basic Flow</b>	1/ System has displayed the chosen photo at the end of <i>View photo</i> use case 2/ User chooses “delete” option to delete the photo 3/ System prompts the user to confirm action 4/ Confirmation of delete action
<b>Exceptions</b>	4/ if System cannot delete the photo 4.1/ System shows an error message 4.2/ Return to step 1
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	
<b>Open Issues:</b>	

<b>Use Case Number:</b>	16
<b>Use Case Name:</b>	Delete Group/Skin Condition
<b>Participating Actors:</b>	main user (patient)
<b>Goal:</b>	Delete some group or some skin condition
<b>Trigger:</b>	User chooses the <i>Delete</i> option in the menu after the System shows a list of groups or a list of skin conditions
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, remove the chosen group or skin condition out of the information of all photos that belong to the group or skin condition
<b>Basic Flow</b>	1/ System has displayed the list of groups or the list of skin conditions at the end of <i>View groups</i> or <i>View skin conditions</i> 2/ User chooses <i>delete</i> option in the menu 3/ System prompts User to choose groups or skin conditions to delete 4/ User chooses groups or skin conditions to delete 5/ System prompts User to confirm deletion 6/ User confirms deletion 7/ System identifies all photos that belong to the deleted group or deleted skin condition and remove the group or skin condition out of the information of those photos; remove the group or skin condition
<b>Exceptions</b>	7/ if System cannot modify info of the photos 7.1/ System shows an error message 7.2/ Not delete the group or skin condition 7.3/ Return to previous activity
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	

<b>Open Issues:</b>	
---------------------	--



<b>Use Case Number:</b>	17
<b>Use Case Name:</b>	View Multiple Photos
<b>Participating Actors:</b>	main User
<b>Goal:</b>	Show the list of photos
<b>Trigger:</b>	The User chooses a criteria to view a list of photos
<b>Precondition:</b>	
<b>Postcondition:</b>	On success, System displays the list of photos or an empty list if there is no photo meet the criteria
<b>Basic Flow</b>	1/ System displays a list of photos meet the criteria 2/ User click on a photo to view it bigger
<b>Exceptions</b>	
<b>Qualities:</b>	
<b>Functionality:</b>	
<b>Constraints:</b>	
<b>Includes:</b>	
<b>Extends:</b>	
<b>Related Artifacts:</b>	
<b>Notes:</b>	
<b>Open Issues:</b>	

-Tests for use case 1-

- 1/ check if the application launches properly
- 2/ check if all the gui and make sure it works
- 3/ check if the application closes with out errors

-Tests for use case 2-

- 1/ verify the application can capture an image with its camera
- 2/ verify the photo can be displayed on screen
- 3/ test to see if all user interaction with the UI works correctly such as user input, and the user selections on the skin condition, and group options.
- 4/ verify that the picture can be stored
- 5/ verify that the picture can be reloaded after being stored
- 6/ verify how the program handles running out of memory to store a new photo

-Tests for use case 3-

- 1/ check that stored photos have a time stamp
- 2/ check that they have the correct time stamp

-Tests for use case 4-

- 1/ verify the user may view a photo

-Tests for use case 5-

- 1/ verify the application can display two photos properly
- 2/ verify the program acts appropriately when the user only has 0 images, and when the user only has 1 image
- 3/ verify the user can choose two photos
- 4/ verify the view works properly for displaying the two photos

-Tests for use case 6-

- 1/ verify that a photo can be added to a group
- 2/ handle the null case, where the user tries to add null to a group?

-Tests for use case 7-

- 1/ verify that the layer works
- 2/ verify a photo can still be capture with the layer
- 3/ handle the case when no photos are present for the user to use as a layer

-Tests for use case 8-

- 1/ verify a reminder can be created

-Tests for use case 9-

- 1/ verify the user is able to use the retake photo option
- 2/ verify that the previous photo has been removed from the application
- 3/ verify the new photo is able to be stored by the application

-Tests for use case 10-

- 1/ verify all the skin conditions are displayed correctly
- 2/ verify the case where no skin conditions exist is handled correctly

-Tests for use case 11-

- 1/ verify all the appropriate groups are displayed
- 2/ verify the case is handled where there are no groups

-Tests for use case 12-

- 1/ verify the case is handled when no skin condition is selected
- 2/ check to see if all the photos with that skin condition are displayed

-Tests for use case 13-

- 1/ verify the case where no groups are selected
- 2/ verify that the correct photos are displayed corresponding to the group

-Tests for use case 14-

- 1/ verify the photo is saved properly after being edited
- 2/ verify the new fields are updated correctly
- 3/ handle the case where the new fields are null

-Tests for use case 15-

- 1/ verify that the photo is removed, and all other data relating to that photo is also removed

-Tests for use case 16-

- 1/ verify that the group/skin condition is deleted
- 2/ verify that all photos who were in the group or had the skin condition tags now no longer have the tags

-Tests for use case 17-

- 1/ verify that the list of photos is displayed correctly
- 2/ verify that the list of photos for 0 photos doesn't create errors